

BD FACSCalibur

Configuration and Fluorophore Guide for FACS Calibur		
Laser, Detector	Filter (nm)	Example Fluorophores (AF= Alexa Fluor)
Blue FL1	530/30	FITC, GFP, CFSE, AF488
Blue FL2	585/42	PE, DSRed
Blue FL3	670 LP	PI, PerCP, PE-Cy5
Red FL4	661/16	APC, AF647

For any enquiries, please contact:

Kylie Price

Flow Cytometry Suite Manager
E: kprice@malaghan.org.nz
P: +64 4 499 6914 ext 850

Alfonso Schmidt

Flow Cytometry Assistant
E: aschmidt@malaghan.org.nz
P: +64 4 499 6914 ext 868