

## FACS Vantage DiVa

Configuration and Fluorophore FACS Vantage DiVa		
Laser, Detector	Filter (nm)	Example Fluorophores (AF= Alexa Fluor)
UV A	450/50	DAPI, Hoechst Blue
UV B	620 LP	Hoechst Red, QDot 800
Blue A	530/30	FITC, AF488, GFP
Blue B	585/42	PE, CTO
Blue C	670 LP	PI, PE-Cy5
Blue D	780/60	PE-Cy7
Red A	661/16	APC, AF647
Red B	780/60	APC-AF750, APC-Cy7
Yellow Green A	594/25	DsRed, dTomato
Yellow Green B	610/20	mCherry

For any enquiries, please contact:

**Kylie Price**

Flow Cytometry Suite Manager  
E: [kprice@malaghan.org.nz](mailto:kprice@malaghan.org.nz)  
P: +64 4 499 6914 ext 850

**Alfonso Schmidt**

Flow Cytometry Assistant  
E: [aschmidt@malaghan.org.nz](mailto:aschmidt@malaghan.org.nz)  
P: +64 4 499 6914 ext 868